


PRODUCT UPDATE

At Fisher & Paykel Healthcare we are continuously reviewing and improving our product range to meet our customers' needs. We are pleased to inform you that the flow rates of Optiflow™ Junior 2 nasal interfaces, when used with F&P 950™ Neonatal circuits, are as follows:







Flow rates shown describe technical capability of the product when used at sea level. Ensure clinical judgement is applied when prescribing flow rates.

Refer to manufacturer's user instructions for correct usage, product setup guidance, warnings, cautions and contraindications.

Please contact your local Fisher & Paykel Healthcare representative for information about this or any of our other products or therapies, or visit www.fphcare.com



950N40 Flow rates (L/min)

	OJR410 XS	0.5–10
	OJR412 S	0.5–10
	OJR414 M	0.5–11
	OJR416 L	0.5–28
	OJR418 XL	0.5–31
	OJR520 XXL	0.5–36